

Precise Diagnosis
at Anytime, Anywhere

SE-310 Electrocardiograph



SE-310

3-Channel Electrocardiograph

Combines lightweight with compact design, SE-310 is the ideal tool for healthcare professionals on the go. Whether in clinics, GP offices, or rehabilitation settings, this portable device offers an instant, precise diagnostics.

In addition, the flat panels feature a smooth surface for effortless cleaning and disinfection. Its minimalist form allows for quick and thorough sanitization, ensuring the highest standards of hygiene are maintained in every scenarios.



Lightweight and Portable Design



Capacitive Touchscreen and short keys



Aluminum handle (Optional)



Extraordinary Performance

CMRR

A high Common Mode Rejection Ratio up to 140 dB restrains signal noises to the greatest extent, thus offers you traces with excellent quality.

Bandwidth

A wide bandwidth from 0.01 to 350 Hz guarantees the accuracy of the ECG signal as well as the interpretation especially for pediatric patients.

Sampling Rate

The leading sampling rate as high as 64,000 Hz aims at capturing all the details to generate precise result.

Signal Quality Detection

By showing different colors to indicate normal, poor or missing signals, this function helps you to identify the signal quality even at a glance.

Reliable ECG Algorithm

SEMIP ECG algorithm offers an accurate and reliable reference to maximize clinicians' confidence in making diagnosis.

On-Screen Diagnosis

Comprehensive diagnostic tools improve clinicians productivity, reduce time spent on diagnosis and let them focus on the patients.

Seamless Connectivity





SE-310

Electrocardiograph



Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-ECG-SE-310
V1.0-20250106