

# H100B

Handheld Pulse Oximeter



Outpatient



Transport



Ward



**EDAN**

A world of potential

# H100B Handheld Pulse Oximeter

As a handheld pulse oximeter, H100B is designed to serve the needs of both long-term monitoring and spot check. It's dedicated to offer reliable, accurate and sensitive measurement of blood oxygen saturation under even harsh conditions.

## Measured Parameters: SpO<sub>2</sub> + PR

- LCD display with backlight control
- Dual work modes: monitoring & spot check
- Numeric display with plethysmogram display
- Trend review
- Pulse-tone modulation (Pitch Tone)
- Adjustable audio and visual alarms
- Powerful data storage capacity (up to 300 hours)
- PatientCare Viewer software for PC data management
- 4 x AA batteries for up to 48 hours of work

**300 h**  
SpO<sub>2</sub> Review

**100**  
Patient ID

- Rechargeable battery for up to 30 hours of work
- Optional battery charger stand
- Automatic power-off function for power saving



## Features



Unique Sensor Design



Bed Rail Mount



Battery Charger Stand



Carrying Bag